

Catalog Number: IRT003

## Anti-DYKDDDDK-tag antibody


Package size: 25, 100 µl  
Store at: -20°C  
MW (kDa): 61

### Overview

<b>Product Name</b>	DYKDDDDK-tag antibody (Binds to same epitope as Sigma's Anti-FLAG® M2 Antibody)
<b>Product Number</b>	IRT003
<b>Clonality</b>	Polyclonal
<b>Host</b>	Rabbit
<b>Immunogen</b>	Synthetic peptide (DYKDDDDK)
<b>Recommended Applications Dilutions</b>	Western Blot 1:3000-1:5000 Immunofluorescence 1:500 – 1:1000
<b>Storage Buffer</b>	100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
<b>Concentration</b>	1.51 mg/ml
<b>Purity</b>	Affinity column purified
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C
<b>Notes</b>	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

### Image

#### IRT003 anti DYKDDDDK-tag antibody WB image



All lanes : Anti DYKDDDDK-tag antibody at 1/3000 dilution

Lysates/proteins at 35 µg per lane

This blot was produced using a 7.5% SDS-PAGE. The gel was run at 140V for 50 minutes before being transferred onto a Nitrocellulose membrane at 18V for 60 minutes. The membrane was then blocked for an hour before being incubated with IRT003 overnight at 4°C.

IReal Biotechnology Co., Ltd.

