

Datasheet for IR30-78


## anti- CD133 / Prominin-1 antibody

Package size: 25, 100  $\mu$ l  
Store at: -20°C  
MW (kDa):

### Overview

<b>Product Name</b>	CD133 / Prominin-1 antibody ( <b>Extracellular domain</b> )
<b>Product Number</b>	IR30-78
<b>Gene Description</b>	prominin 1
<b>Clonality</b>	Polyclonal
<b>Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Recommended Applications Dilutions</b>	Western Blot <b>Assay-dependent</b> Immunofluorescence 1:200 Immunohistochemistry (Paraffin) 1:100
<b>Storage Buffer</b>	100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
<b>Concentration</b>	0.71 mg/ml
<b>Purity</b>	Affinity column purified
<b>Storage</b>	Store at +4°C for short term storage. Long time storage is recommended at -20°C
<b>Notes</b>	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

### IR30-78 anti- CD133 / Prominin-1 antibody IHC image



Immunohistochemical analysis of paraffin embedded Human cancer tissue labeling CD133 / Prominin-1 with IR30-78 at 1/100.

IReal Biotechnology Co., Ltd.

7F., No. 9, Sec. 3, Zhonghua Rd., North Dist., Hsinchu City 300, Taiwan (R.O.C.)

Website : [www.irealbio.com](http://www.irealbio.com)

TEL : 03-526-0005

E-mail : [irealbio@irealbio.com](mailto:irealbio@irealbio.com)