

Catalog Number: IR254-903

Anti-CD16 antibody


Package size: 25, 100 μ l
Store at: -20°C
MW (kDa): -

Overview

Product Name	CD16 antibody
Product Number	IR254-903
Gene Description	CD16 molecule
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human
Recommended Applications Dilutions	Flow cytometry 1:100 – 1:200
Storage Buffer	100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1.62 mg/ml
Purity	Affinity column purified
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

Image

IR254-903 anti-CD16 antibody FACS image



Product: CD16 antibody
Catalog: IR254-903
Cell: Human peripheral blood lymphocytes
CD16 antibody: 1-100
Data from: POLOWANG BIOTECH INC

IReal Biotechnology Co., Ltd.

