

Catalog Number: IR252-897-FITC


## Anti-CD3 (FITC conjugated) antibody

Package size: 25, 100  $\mu$ l  
 Store at: 4°C  
 MW (kDa): -

### Overview

Product Name	CD3 (FITC conjugated) antibody
Product Number	IR252-897-FITC
Gene Description	CD3e molecule, epsilon (CD3-TCR complex)
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human
Recommended Applications Dilutions	Immunofluorescence 1:200 Flow Cytometry 1:100 – 1:200
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN <sub>3</sub> was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored undiluted between 2°C and 8°C, protected from prolonged exposure to light. <b>Do not freeze.</b>
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

### IR252-897-FITC anti-CD3 (FITC conjugated) antibody FACS image



Product: CD3 (FITC conjugated) antibody  
 Catalog: IR252-897-FITC  
 Cell: Human peripheral blood  
 CD3 antibody: 1-100  
 Data from: POLOWANG BIOTECH INC

iReal Biotechnology Co., Ltd.

