

Catalog Number: IR208-752-FITC


Anti-CD44 (FITC conjugated) antibody

Package size: 25, 100 μ l
 Store at: 4°C
 MW (kDa): -

Overview

Product Name	CD44 (FITC conjugated) antibody
Product Number	IR208-752-FITC
Gene Description	CD44 molecule
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human, mouse, rat
Recommended Applications Dilutions	Immunofluorescence 1:200 Flow Cytometry 1:100 – 1:200
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN ₃ was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored undiluted between 2°C and 8°C, protected from prolonged exposure to light. Do not freeze.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IR208-752-FITC anti-CD44 (FITC conjugated) antibody FACS image



Product: CD44 (FITC conjugated) antibody
 Catalog: IR208-752-FITC
 Cell: MDA-MB-231
 CD44 antibody: 1-100
 Data from: POLOWANG BIOTECH INC

IReal Biotechnology Co., Ltd.

