

Catalog Number: IR208-752-F647


## Anti-CD44 (**HiLyte Fluor647 conjugated**) antibody

Package size: 25, 100  $\mu$ l  
Store at: 4°C  
MW (kDa): -

### Overview

Product Name	CD44 ( <b>HiLyte Fluor647 conjugated</b> ) antibody
Product Number	IR208-752-F647
Gene Description	CD44 molecule
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human, mouse, rat
Recommended Applications Dilutions	Immunofluorescence 1:200 Flow Cytometry 1:100 – 1:200
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN <sub>3</sub> was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored undiluted between 2°C and 8°C, protected from prolonged exposure to light. <b>Do not freeze.</b>
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

### IR208-752-F647 anti-CD44 (**HiLyte Fluor647 conjugated**) antibody FACS image



Product: CD44 (**HiLyte Fluor647 conjugated**) antibody  
Catalog: IR208-752-F647  
Cell: MDA-MB-231  
CD44 antibody: 1-150  
Data from: POLOWANG BIOTECH INC

iReal Biotechnology Co., Ltd.

