

Catalog Number: IR20-56-FITC


Anti-NeuN (FITC conjugated) antibody

Package size: 25, 100 μ l
 Store at: 4°C
 MW (kDa): -

Overview

Product Name	NeuN (FITC conjugated) antibody
Product Number	IR20-56-FITC
Gene Description	RNA binding protein, fox-1 homolog (C. elegans) 3
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human, mouse, rat
Recommended Applications Dilutions	Immunofluorescence 1:200 Immunohistochemistry (Frozen) 1:100 – 1:200
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN ₃ was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored undiluted between 2°C and 8°C, protected from prolonged exposure to light. Do not freeze.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IR20-56-FITC anti-NeuN (FITC conjugated) antibody IF image



Immunohistochemical analysis.
 Sample: Frozen-sectioned mouse cerebellum.
 NeuN protein stained by NeuN antibody (IR20-56-FITC conjugated) diluted at 1:100.

IReal Biotechnology Co., Ltd.