

Catalog Number: IR218-3-F647

Anti- EpCAM extracellular (**HiLyte Fluor647 conjugated**) Antibody


Package size: 100 μ l
Store at: 4°C

Overview

Product Name	EpCAM extracellular (HiLyte Fluor647 conjugated) antibody
Product Number	IR218-3-F647
Gene Description	epithelial cell adhesion molecule
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human
Recommended Applications Dilutions	Immunofluorescence 1:100 – 1:200 Flow Cytometry 1:100
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN3 was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored at 4°C. Do not aliquot the antibody. Protect from light. Do not freeze.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

Image

IR218-3-F647 anti-EpCAM extracellular (**HiLyte Fluor647 conjugated**) antibody FACS image



Product: EpCAM extracellular (**HiLyte Fluor647 conjugated**) antibody
Catalog: IR218-3-F647
Cell: MCF-7
EpCAM extracellular antibody: 1-100
Data from: POLOWANG BIOTECH INC

IReal Biotechnology Co., Ltd.