

Catalog Number: IR329-1034

Anti Toll-like Receptor 9 / TLR9 antibody

Package size: 25, 100 μ l


Store at: -20°C

MW (kDa): 130


Overview

Product Name	Toll-like Receptor 9 / TLR9 antibody
Product Number	IR329-1034
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Recommended Applications Dilutions	Western Blot 1:500- 1:1000 Immunofluorescence 1:200– 1:400 Immunohistochemistry (Paraffin) 1:100 – 1:150
Storage Buffer	100mM Tris Glycine, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1.5 mg/ml
Purity	Affinity column purified
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IR329-1034 Toll-like Receptor 9 / TLR9 antibody WB image



IR329-1034 Toll-like Receptor 9 / TLR9 antibody IHC image



Toll-like Receptor 9 / TLR9 antibody - at 1/1000 dilution
Lysates at 40 μ g per lane. This blot was produced using a 5% SDS-PAGE. NC membrane was then blocked for an hour before being incubated with IR329-1034 overnight at 4°C.

Immunohistochemical analysis of paraffin embedded Human tonsil cancer tissue labeling **Toll-like Receptor 9 / TLR9** with IR329-1034 at 1/150.