

Catalog Number: IRPA2301


Anti-Cyclophilin A / PPIA antibody

Package size: 25, 100 μ l
Store at: -20°C
MW (kDa): 17

Overview

Product Name	Cyclophilin A / PPIA antibody
Product Number	IRPA2301
Gene Description	peptidylprolyl isomerase A (cyclophilin A)
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Recommended Applications Dilutions	Western Blot 1:1000 – 1:2000 Immunofluorescence 1:300 – 1:500
Storage Buffer	100mM Tris Glycine, 35% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRPA2301 Cyclophilin A / PPIA antibody WB image



IReal Biotechnology Co., Ltd.