

Catalog Number: IRV008

Anti-Norovirus VP1 Monoclonal antibody [IR11C]


Package size: 25, 100 μ l
Store at: -20°C
MW (kDa):

Overview

Product Name	Norovirus VP1 Monoclonal antibody [IR11C]
Product Number	IRV008
Clonality	Monoclonal
Host	Mouse
Antigen Species	Virus
Recommended Applications Dilutions	Western Blot 1:3000-1:15000 Immunofluorescence 1:500 – 1:1000 ELISA (1:5000-10000)
Storage Buffer	100mM Tris Glycine, 30% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRV008 anti Norovirus VP1 Monoclonal antibody [IR11C] WB image

Norovirus VP1
Recombinant protein



Anti- Norovirus VP1 Monoclonal antibody [3-8A11C] - at 1/3000 dilution

Recombinant proteins at 0.1 μ g per lane.
Nitrocellulose membrane was then blocked for an hour before being incubated with IRV008 overnight at 4°C.

iReal Biotechnology Co., Ltd.