

Catalog Number: IRM323

Anti IL-2 antibody [MB2]


Package size: 25, 100 μ l
Store at: -20°C
MW (kDa): -

Overview

Product Name	IL-2 antibody [MB2]
Product Number	IRM323
Clonality	Monoclonal
Host	Mouse (mouse IgG1)
Species Reactivity	Human
Recommended Applications Dilutions	Western Blot 1:500 Immunofluorescence 1:100 – 1:300
Storage Buffer	Phosphate-buffered solutions, pH 7.2, 50% Glycerol, 0.025% ProClin 300 was added as a preservative
Concentration	1mg/ml
Purity	Protein A column purified
Storage	Long time storage is recommended at -20°C
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

Image

IRM323 IL-2 antibody [MB2] WB Image



Anti IL-2 antibody [MB2] at 1/500 dilution

Recombinant Human Interleukin-2 (hIL-2) proteins. Nitrocellulose membrane was then blocked for an hour before being incubated with IRM323 overnight at 4°C.

