

Catalog Number: IRR006


Anti-mouse CD8a antibody [53-6.7]

Package size: 100 μ l
Store at: -20°C
MW (kDa): -

Overview


Product Name	mouse CD8a antibody [53-6.7]
Product Number	IRR006
Clonality	Monoclonal
Host	Rat (Rat IgG2a)
Species Reactivity	Mouse
Recommended Applications Dilutions	Immunohistochemistry 1:100 Flow Cytometry 1:50 – 1:100
Storage Buffer	Phosphate-buffered solution, pH 7.2, containing 0.025% ProClin 300, 20% Glycerol
Concentration	0.5 mg/ml
Purity	Affinity column purified
Storage	Store at -20°C.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRR006 mouse CD8a antibody [53-6.7] IHC image




Immunohistochemical analysis of paraffin embedded mouse spleen tissue labeling mouse CD8a antibody [53-6.7] with IRR006 at 1/100.

IRR006 mouse CD8a antibody [53-6.7] FACS image



Mouse peripheral blood monocytes stained with IRR006 mouse CD8a antibody [53-6.7] antibody.

IRR006 mouse CD8a antibody [53-6.7] IHC-Fr image



CD8a was detected in immersion fixed mouse spleen using Rat Anti-Mouse CD8a antibody [53-6.7] (IRR006) at 1-100, at room temperature. counterstained with DAPI (blue).

