

Catalog Number: IRR007

Anti-F4/80 [BM8] monoclonal antibody


Package size: 100 μ l
Store at: -20°C
MW (kDa): -

Overview


Product Name	F4/80 [BM8] monoclonal antibody
Product Number	IRR007
Clonality	Monoclonal
Host	Rat (Rat IgG2a)
Species Reactivity	Mouse
Recommended Applications Dilutions	Immunohistochemistry 1:100 Flow Cytometry 1:50 – 1:100
Storage Buffer	Phosphate-buffered solution, pH 7.2, containing 0.025% ProClin 300, 20% Glycerol
Concentration	0.5 mg/ml
Purity	Affinity column purified
Storage	Store at -20°C.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IRR007 F4/80 [BM8] monoclonal IHC image

IRR007 F4/80 [BM8] monoclonal FACS image



Mouse brain



Immunohistochemical analysis of paraffin embedded mouse brain tissue labeling F4/80 [BM8] with IRR007 at 1/100.

Mouse peritoneal macrophages stained with F4/80 [BM8] antibody (IRR007), followed by anti-Rat-FITC

IReal Biotechnology Co., Ltd.

