

Influenza A Virus Nucleoprotein

DATA SHEET

Product Name: Influenza A Virus Nucleoprotein

Cat. No: IRRP7010

Product type: Protein


Description: Influenza A Virus Nucleoprotein produced in *Escherichia coli*

Source: *Escherichia Coli*.

Solubility: Reconstitution with PBS.

Purity: >95% as determined by SDS-PAGE

Stability&Storage: Samples are stable for up to twelve months from date of receipt at -20°C to -80°C. Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.



Notes: *For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.*