

## Influenza B Virus Nucleoprotein

## DATA SHEET

**Product Name:** Influenza B Virus Nucleoprotein

**Cat. No:** IRRP7011

**Product type:** Protein


**Description:** Influenza B Virus Nucleoprotein produced in Escherichia coli

**Source:** *Escherichia Coli*.

**Solubility:** Reconstitution with PBS.

**Purity:** >95% as determined by SDS-PAGE

**Stability&Storage:** Samples are stable for up to twelve months from date of receipt at -20°C to -80°C. Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.



**Notes:** *For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.*