

Catalog Number: IRFS4132


**Mouse anti-rabbit IgG, Texas red conjugated**

Package size: 200, 100  $\mu$ l  
Store at: 4°C

### Overview

|                                    |  |
|------------------------------------|--|
| Product Name                       | Mouse anti-rabbit IgG, Texas red conjugated  |
| Product Number                     | IRFS4132   |
| Label                              | Texas red conjugated   |
| Clonality                          | Monoclonal   |
| Host                               | Mouse  |
| Species Reactivity                 | Rabbit IgG   |
| Recommended Applications Dilutions | Immunofluorescence 1:200 – 1:400<br>flow cytometry (0.5-1 $\mu$ g per 1 x 10 <sup>6</sup> cells)<br>Optimal dilution to be determined by titration.<br>(Other Assays User Optimized) |
| Storage Buffer                     | PBS containing 1% stabilizer protein and 0.025% ProClin 300.   |
| Concentration                      | 1 mg/ml  |
| Storage                            | Store at 4° C, <b>**DO NOT FREEZE**</b> . Stable for one year from the date of shipment. Non-hazardous. No MSDS required.  |
| SOURCE                             | Mouse anti-rabbit IgG-Texas red is an affinity purified secondary antibody raised in mouse against rabbit IgG and conjugated to Texas red.   |

### Image



**IRFS4132. Mouse anti-rabbit IgG, Texas red conjugated ICC/IF image.**

Primary antibody : Histone H3 antibody, Catalog : IR4-10.  
Cell: A431.

Detection reagent used: Mouse anti-rabbit IgG, Texas red conjugated (IRFS4132).

