

Catalog Number: IRFS4142


**Mouse anti-rabbit IgG, APC conjugated**

Package size: 200, 100  $\mu$ l  
Store at: 4°C

## Overview

Product Name	Mouse anti-rabbit IgG, APC conjugated
Product Number	IRFS4142
Label	APC conjugated
Clonality	Monoclonal
Host	Mouse
Species Reactivity	Rabbit IgG
Recommended Applications Dilutions	Immunofluorescence 1:200 – 1:400 flow cytometry (0.5-1 $\mu$ g per 1 x 10 <sup>6</sup> cells) Optimal dilution to be determined by titration. (Other Assays User Optimized)
Storage Buffer	PBS containing 1% stabilizer protein and 0.025% ProClin 300.
Concentration	1 mg/ml
Storage	Store at 4° C, <b>**DO NOT FREEZE**</b> . Stable for one year from the date of shipment. Non-hazardous. No MSDS required.
SOURCE	Mouse anti-rabbit IgG-APC is an affinity purified secondary antibody raised in mouse against rabbit IgG and conjugated to APC.

## Image



**IRFS4142. Mouse anti-rabbit IgG, APC conjugated ICC/IF image.**

Primary antibody : Histone H3 antibody, Catalog : IR4-10.  
Cell: A431.

Detection reagent used: Mouse anti-rabbit IgG, APC conjugated (IRFS4142).

