

Catalog Number: IR75-253-F647


Anti-mouse CD31 (HiLyte Fluor647 conjugated) antibody

Package size: 25, 100 μ l
Store at: 4°C
MW (kDa): -

Overview

Product Name	Mouse CD31 (HiLyte Fluor647 conjugated) antibody
Product Number	IR75-253-F647
Clonality	Polyclonal
Host	Rabbit
Species Reactivity	Mouse specific
Recommended Applications Dilutions	Immunofluorescence 1:200 – 1:500 Flow Cytometry 1:100 – 1:200
Storage Buffer	PBS (pH 7.4), 0.2% BSA, 0.1% Gelatin, 0.1% NaN ₃ was added as a preservative
Concentration	1 mg/ml
Purity	Affinity column purified
Storage	This product is stable for 12 months when stored undiluted between 2°C and 8°C, protected from prolonged exposure to light. Do not freeze.
Notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

IR75-253-F647 anti-mouse CD31 (HiLyte Fluor647 conjugated) antibody FACS image



Product: mouse CD31 (HiLyte Fluor647 conjugated) antibody.
Catalog: IR75-253-F647
Sample: C57BL/6 mouse splenocytes
Mouse CD31 antibody: 1-150
Data from: POLOWANG BIOTECH INC

IReal Biotechnology Co., Ltd.

