

Neuronal Specific Enolase Recombinant Protein Data Sheet

Cat No: IRRP7006

Product Name: NSE/ENO2 recombinant protein

Protein Construction: A DNA sequence encoding the mature form of human ENO2 (P09104)(Met 1-Leu 434) was fused with a polyhistidine tag at the N-terminus.

Concentrate: 0.1 mg/ml

Vector: pET28

Expressed Host: *E.coli*

Form Supplied: Liquid (100µg)

Storage buffer: 100 mM tris pH7.4, 200 mM NaCl,

Storage instructions: Aliquot and store at -20°C.

